

### 383: DESULFOBACTERIUM MEDIUM

<b>Solution A</b>	952.00	ml
<b>Solution B</b>	20.00	ml
<b>Solution C</b>	20.00	ml
<b>Solution D</b>	1.00	ml
<b>Solution E</b>	10.00	ml

1. Solution A is sparged with 80% N<sub>2</sub> and 20% CO<sub>2</sub> gas mixture to reach a pH below 6 (at least 30 min), then distributed in anoxic cultivation vials and autoclaved under the same gas atmosphere. Solution B is autoclaved separately under 80% N<sub>2</sub> and 20% CO<sub>2</sub> gas atmosphere. Solutions C and D are prepared under 100% N<sub>2</sub> gas and filter-sterilized. Solution E is autoclaved under 100% N<sub>2</sub> gas. To complete the medium appropriate amounts of solutions B to E are added to the sterile solution A in the sequence as indicated. Final pH of the medium should be 7.0 – 7.2.

2. Note: Addition of 10 - 20 mg sodium dithionite per liter (e.g. from 5% (w/v) solution freshly prepared under N<sub>2</sub> and filter-sterilized) may stimulate growth of some strains at the beginning. For transfers use 5 - 10% inoculum. Incubate all strains in the dark.

For DSM 2056: Replace pyruvate with 0.70 g/l Na-butyrate, 0.30 g/l Na-caproate, and 0.15 g/l sodium octanoate.

For DSM 3383: Replace pyruvate with indole as substrate. Dissolve 0.30 g of indole in 90 ml water by heating and shaking in a closed bottle under 100% N<sub>2</sub> gas atmosphere, autoclave, then add 7.00 ml of a sterile anoxic stock solution of NaCl (30% w/v) and 0.70 ml of a sterile anoxic stock solution of MgCl<sub>2</sub> × 6 H<sub>2</sub>O (40% w/v). Store the indole-salt solution in the dark. Reheat and shake before use. Add to sterile medium 30.00 ml/l of the indole-salt solution in the beginning, and 2 x 30.00 ml/l during growth.

For DSM 3384: The vitamin solution of medium 120 is replaced with 1.00 ml/l of seven vitamins solution of medium 503. Replace pyruvate with 0.40 g/l Na-benzoate and add 10.00 ml/l clarified rumen fluid (see medium 1310) as growth supplement. Initiation of growth requires addition of freshly prepared sodium dithionite.

For DSM 4661: Replace pyruvate with resorcinol. During growth the culture is fed once with the same amount of resorcinol.

For DSM 5091, DSM 9788: Replace pyruvate with 2.00 g/l malonic acid.

For DSM 7044, DSM 7120, DSM 7467, DSM 12567, DSM 13228, DSM 23484: Supplement medium with 1.00 ml/l seven vitamins solution (see medium 503). Replace pyruvate with 0.40 g/l benzoate and 0.10 g/l yeast extract sterilized separately by filtration and add to the autoclaved medium from anoxic stock solutions.

For DSM 8540: Replace pyruvate with 0.30 g/l p-hydroxybenzoate.

For DSM 9705: Replace pyruvate with 1.00 g/l Na-glycolate sterilized by filtration.

## 383: DESULFOBACTERIUM MEDIUM

For DSM 12861, DSM 12883: Replace pyruvate with 1.00 g/l valeric acid.

For DSM 12888: Replace pyruvate with 1.10 g/l Na-butyrate and 1.10 g/l caproic acid.

For DSM 13036: Replace pyruvate with 1.00 g/l betaine.

For DSM 14454: Naphthalene dissolved in heptamethylnonane can be used as substrate instead of pyruvate.

For DSM 15576, DSM 16219: Replace pyruvate with 1.00 g/l Na-caprylate and supplement medium with 1.00 ml/l seven vitamins solution (see medium 503) added to the autoclaved medium from an anoxic stock solution sterilized by filtration.

For DSM 17291: Replace pyruvate with 3.00 g/l Casamino acids.

For DSM 17477: Replace pyruvate with 1.90 g/l sodium glutamate monohydrate and 1.00 g/l yeast extract.

For DSM 21856: Replace pyruvate with 0.90 ml/l 1-butanol.

For DSM 28890: Replace pyruvate with 0.50 g/l yeast extract and 0.68 g/l sodium formate.

For DSM 100305, DSM 105015: Change amount of pyruvate to 2.20 g/l and add 0.10 g/l yeast extract.

### Solution A

Na <sub>2</sub> SO <sub>4</sub>	3.00	g
KH <sub>2</sub> PO <sub>4</sub>	0.20	g
NH <sub>4</sub> Cl	0.30	g
NaCl	21.00	g
MgCl <sub>2</sub> x 6 H <sub>2</sub> O	3.00	g
KCl	0.50	g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	0.15	g
<b>Trace element solution SL-10</b>	1.00	ml
<b>Selenite-tungstate solution</b>	1.00	ml
Sodium resazurin (0.1% w/v)	0.50	ml
Distilled water	950.00	ml

### Solution B

Na <sub>2</sub> CO <sub>3</sub>	1.00	g
Distilled water	20.00	ml

### Solution C

Na-pyruvate	2.50	g
Distilled water	20.00	ml

### Solution D

## 383: DESULFOBACTERIUM MEDIUM

<b>Wolin's vitamin solution (10x)</b>	1.00	ml
---------------------------------------	------	----

### Solution E

Na <sub>2</sub> S x 9 H <sub>2</sub> O	0.40	g
Distilled water	10.00	ml

### Selenite-tungstate solution (from medium 385)

NaOH	0.50	g
Na <sub>2</sub> SeO <sub>3</sub> x 5 H <sub>2</sub> O	3.00	mg
Na <sub>2</sub> WO <sub>4</sub> x 2 H <sub>2</sub> O	4.00	mg
Distilled water	1000.00	ml

### Trace element solution SL-10 (from medium 320)

HCl (25%)	10.00	ml
FeCl <sub>2</sub> x 4 H <sub>2</sub> O	1.50	g
ZnCl <sub>2</sub>	70.00	mg
MnCl <sub>2</sub> x 4 H <sub>2</sub> O	100.00	mg
H <sub>3</sub> BO <sub>3</sub>	6.00	mg
CoCl <sub>2</sub> x 6 H <sub>2</sub> O	190.00	mg
CuCl <sub>2</sub> x 2 H <sub>2</sub> O	2.00	mg
NiCl <sub>2</sub> x 6 H <sub>2</sub> O	24.00	mg
Na <sub>2</sub> MoO <sub>4</sub> x 2 H <sub>2</sub> O	36.00	mg
Distilled water	990.00	ml

First dissolve FeCl<sub>2</sub> in the HCl, then dilute in water, add and dissolve the other salts. Finally make up to 1000.00 ml.

### Wolin's vitamin solution (10x) (from medium 120)

Biotin	20.00	mg
Folic acid	20.00	mg
Pyridoxine hydrochloride	100.00	mg
Thiamine HCl	50.00	mg
Riboflavin	50.00	mg
Nicotinic acid	50.00	mg
Calcium D-(+)-pantothenate	50.00	mg
Vitamin B <sub>12</sub>	1.00	mg
p-Aminobenzoic acid	50.00	mg
(DL)-alpha-Lipoic acid	50.00	mg
Distilled water	1000.00	ml