

**520. RUMINOCLOSTRIDIUM CELLULOLYTICUM (CM3) MEDIUM**

(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	1.30	g
KH <sub>2</sub> PO <sub>4</sub>	1.50	g
K <sub>2</sub> HPO <sub>4</sub> x 3 H <sub>2</sub> O	2.90	g
FeSO <sub>4</sub> x 7 H <sub>2</sub> O solution (0.1% w/v in 0.1 N H <sub>2</sub> SO <sub>4</sub> )	1.25	ml
Trace element solution SL-10 (see medium 320)	1.00	ml
Yeast extract	2.00	g
Na-resazurin solution (0.1% w/v)	0.50	ml
MgCl <sub>2</sub> x 6 H <sub>2</sub> O	0.20	g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	75.00	mg
Cellobiose	6.00	g
L-Cysteine-HCl x H <sub>2</sub> O	0.50	g
Na <sub>2</sub> CO <sub>3</sub>	2.50	g
Distilled water	1000.00	ml

Dissolve ingredients except magnesium chloride, calcium chloride, cellobiose, cysteine and carbonate, then sparge medium with 80% N<sub>2</sub> and 20% CO<sub>2</sub> gas mixture for 30 – 45 min to make it anoxic. Dispense medium under the same gas atmosphere into anoxic Hungate-type tubes or serum vials and autoclave. After autoclaving add magnesium chloride, calcium chloride and cellobiose from anoxic stock solutions prepared under 100% N<sub>2</sub> gas. Cellobiose has to be sterilized by filtration. Prior to inoculation add cysteine from a sterile anoxic stock solution prepared under 100% N<sub>2</sub> gas and adjust pH to 7.2 by adding a sterile anoxic stock solution of Na<sub>2</sub>CO<sub>3</sub> (5% w/v) prepared under 80% N<sub>2</sub> and 20% CO<sub>2</sub> gas atmosphere.

*Note: Some strains can be adapted to cellulose as substrate using 10.00 g/l cellulose powder MN 301 (MACHEREY-NAGEL).*

For [DSM 5974](#), [DSM 9801](#) and [DSM 17427](#) use 5.00 g/l D-glucose as the substrate. Supplement medium after autoclaving with 0.50 g/l Na<sub>2</sub>S x 9 H<sub>2</sub>O added from a sterile anoxic stock solution prepared under 100% N<sub>2</sub> gas and adjust pH of complete medium to 6.8 - 7.0, if necessary.