

## 825: METHANOBACTERIUM II MEDIUM

Final pH: 6.8 - 7.0

Final volume: 1011 ml

NH <sub>4</sub> Cl	1.00	g
K <sub>2</sub> HPO <sub>4</sub>	0.30	g
KH <sub>2</sub> PO <sub>4</sub>	0.30	g
MgCl <sub>2</sub> x 6 H <sub>2</sub> O	0.20	g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	0.10	g
KCl	0.10	g
NaCl	0.60	g
<b>Modified Wolin's mineral solution</b>	10.00	ml
Na-acetate	0.50	g
Yeast extract	1.00	g
Sodium resazurin (0.1% w/v)	0.50	ml
NaHCO <sub>3</sub>	4.00	g
<b>Wolin's vitamin solution (10x)</b>	1.00	ml
L-Cysteine HCl x H <sub>2</sub> O	0.50	g
Na <sub>2</sub> S x 9 H <sub>2</sub> O	0.50	g
Distilled water	1000.00	ml

1. Dissolve ingredients except bicarbonate, vitamins, cysteine and sulfide. Sparge medium with 80% H<sub>2</sub> and 20% CO<sub>2</sub> gas mixture for 30 - 45 min to make it anoxic. Add and dissolve bicarbonate, adjust pH to 7.0, then dispense medium under 80% H<sub>2</sub> and 20% CO<sub>2</sub> gas atmosphere into anoxic Hungate-type tubes or serum vials to 30% of their volume and autoclave. Add vitamins (sterilized by filtration), cysteine and sulfide from sterile anoxic stock solutions prepared under 100% N<sub>2</sub> gas. Prior to use check pH of complete medium and adjust to 6.8 - 7.0, if necessary.

2. After inoculation pressurize cultivation vessels to 2 bar overpressure with sterile 80% H<sub>2</sub> and 20% CO<sub>2</sub> gas mixture.

For DSM 24406: Supplement medium with 0.5 g/l DL-dithiothreitol (DTT), 20 ml/l methanol (50% v/v), and 0.5 g/l 2-mercaptopethanesulfonate (coenzyme M) added to the autoclaved medium from sterile anoxic stock solutions. The stock solution of DTT has to be sterilized by filtration. Omit Na<sub>2</sub>S x 9 H<sub>2</sub>O and L-cysteine HCl x H<sub>2</sub>O. Adjust the pH of the final medium to 6.5.

### Modified Wolin's mineral solution (from medium 141)

Nitrolotriacetic acid	1.50	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	3.00	g
MnSO <sub>4</sub> x H <sub>2</sub> O	0.50	g

NaCl	1.00	g
FeSO <sub>4</sub> x 7 H <sub>2</sub> O	0.10	g
CoSO <sub>4</sub> x 7 H <sub>2</sub> O	0.18	g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	0.10	g
ZnSO <sub>4</sub> x 7 H <sub>2</sub> O	0.18	g
CuSO <sub>4</sub> x 5 H <sub>2</sub> O	0.01	g
AlK(SO <sub>4</sub> ) <sub>2</sub> x 12 H <sub>2</sub> O	0.02	g
H <sub>3</sub> BO <sub>3</sub>	0.01	g
Na <sub>2</sub> MoO <sub>4</sub> x 2 H <sub>2</sub> O	0.01	g
NiCl <sub>2</sub> x 6 H <sub>2</sub> O	0.03	g
Na <sub>2</sub> SeO <sub>3</sub> x 5 H <sub>2</sub> O	0.30	mg
Na <sub>2</sub> WO <sub>4</sub> x 2 H <sub>2</sub> O	0.40	mg
Distilled water	1000.00	ml

First dissolve nitrilotriacetic acid and adjust pH to 6.5 with KOH, then add minerals. Adjust final to pH 7.0 with KOH.

### Wolin's vitamin solution (10x) (from medium 120)

Biotin	20.00	mg
Folic acid	20.00	mg
Pyridoxine hydrochloride	100.00	mg
Thiamine HCl	50.00	mg
Riboflavin	50.00	mg
Nicotinic acid	50.00	mg
Calcium D-(+)-pantothenate	50.00	mg
Vitamin B <sub>12</sub>	1.00	mg
p-Aminobenzoic acid	50.00	mg
(DL)-alpha-Lipoic acid	50.00	mg
Distilled water	1000.00	ml