

## Tomato yellow leaf curl virus (TYLCV)

This ELISA to detect TYLCV is a triple antibody sandwich format, TAS-ELISA. It employs a polyclonal antibody (IgG) for coating and a monoclonal antibody (MAb) (developed as broad spectrum reagent against all isolates of TYLCV and TYLCSV) for decorating the virus coat protein. Since the monoclonal antibody is not labelled, a secondary, animal species (mouse) specific antibody is used to react with the bound MAb. This anti mouse (RAM) antibody is labelled with alkaline phosphatase (AP) as reporter group.

Cross reactions with other begomoviruses possible!

Our ELISA reagents are optimized using greiner bio-one microplates, medium binding. Before opening the tubes containing coating antibody (IgG), MAb and RAM-AP- Conjugate please spin down all the liquid by a short centrifugation (approx. 3000rpm for a few seconds).



1. Dilute specific antibody in coating buffer (recommended dilution see delivery note and tube); i.e. 20µl in 20 ml buffer at a recommended dilution of 1:1000 or 40µl in 20 ml buffer at a recommended dilution of 1:500. Add 200µl to each well of the microtiter plate.
2. Cover the plate and incubate at 37 °C for 2- 4 h.
3. Wash plate with PBS-Tween using wash bottle, soak for a few minutes and repeat washing two times. Blot plates by tapping upside down on tissue paper.
4. Add 200 µl of 2% skim milk in PBS-Tween to each well (**blocking**). Cover the plate and incubate for 30 min at 37°C




5. Remove blocking solution and tap dry.
6. Extract samples 1/20 (w/v) in Begomovirus sample extraction buffer (0.05 M Tris containing 0.06 M sodium sulfite, pH 8.5) by grinding young leave tips in extraction buffer.



7. Add 200 µl aliquots of the test sample to duplicate wells.
8. Cover the plate and incubate overnight at 4 °C.
9. Wash three times as in step 3.
10. Add 200 µl of the MAb in appropriate (recommended dilution see delivery note and tube label) conjugate buffer to each well.
11. Cover the plate and incubate at 37 °C for 2- 4 hours.
12. Wash three times as in step 3.



13. Dilute RAM-AP 1:1000 in conjugate buffer, i.e., 20 µl in 20 ml buffer. Add 200 µl to each well.
14. Cover the plate and incubate at 37 °C for 1 hour.
15. Wash three times as in step 3.



17. Add 200 µl aliquots of freshly prepared substrate (1 mg/ ml para-nitrophenyl- phosphate in substrate buffer) to each well.
18. Cover the plate and incubate at 37°C for 30-60 min, or as long as necessary to obtain clear reactions
19. Assess results by:
  - a) Visual observation
  - b) Spectrophotometric measurement of absorbance at 405 nm

#### Reference

Clark, M. F. and Adams. A. N. 1977. Characteristics of the microplate method of enzyme-linked immunosorbent assay for the detection of plant viruses. *Journal of General Virology* 34: 475-483.

### Buffers used in ELISA

#### 1. Coating buffer (pH 9.6)

1.59 g sodium carbonate ( $\text{Na}_2\text{CO}_3$ )  
 2.93 g sodium bicarbonate ( $\text{NaHCO}_3$ )  
 0.20 g sodium azide ( $\text{NaN}_3$ )  
 Dissolve in 900 ml  $\text{H}_2\text{O}$ , adjust pH to 9.6 with HCl and make up to 1 l.

#### 2. PBS (pH 7.4) phosphate buffered saline

8.0 g sodium chloride ( $\text{NaCl}$ )  
 0.2 g monobasic potassium phosphate ( $\text{KH}_2\text{PO}_4$ )  
 1.15 g dibasic sodium phosphate ( $\text{Na}_2\text{HPO}_4$ )  
 0.2 g potassium chloride ( $\text{KCl}$ )  
 0.2 g sodium azide ( $\text{NaN}_3$ )  
 Dissolve in 900 ml  $\text{H}_2\text{O}$ , adjust pH to 7.4 with NaOH or HCl and make up to 1 l.

#### 3. PBS-Tween (PBST)

PBS + 0.5 ml Tween 20 per liter

#### 4. Sample extraction buffer (pH 8.5) for Begomoviruses

0.05 M Tris containing 0.06 M sodium sulfite, pH 8.5

## 5. *Conjugate buffer*

PBST + 2% PVP + 0.2% egg albumin (e.g. Sigma A-5253)

## 6. *Substrate buffer*

97 ml diethanolamine

600 ml H<sub>2</sub>O

0.2 g sodium azide (NaN<sub>3</sub>)

*Adjust to pH 9.8 with HCl and make up to 1 liter with H<sub>2</sub>O*

**Buffers can be stored at 4 ° C for at least 2 months. Warm to room temperature before use.**

## ELISA Troubleshooting

### 1. *No color development*

- a) Did you omit any steps?
- b) Did you use the correct buffer for each step?
- c) Is your enzyme OK? Serum OK?
- d) Is your positive control homologous to antiserum (IgG)?

**Recommendations** - Do a titration plate. Use reliable positive control in each plate. Pretest enzyme conjugate on substrate. Check purified IgG for antigen-specific IgG.

### 2. *Nonspecific color development*

- a) If in edge wells only:
  - Make sure the humidity in the incubator is sufficiently high.
  - If this does not help, don't use edge or border wells, fill with buffer only.
- b) If in whole plate:
  - incomplete washing
  - old substrate
  - use recommended ELISA plate (greiner medium binding)
  - error in loading sequence

**Recommendations** - Use reliable negative control in each plate. Use fresh substrate and check for spontaneous color change. Cover plates while incubating. Check pH of the buffers used.

- c) Some wells with inconsistent or unexpected reactions
  - incomplete washing
  - error in loading test antigens
  - spillage between wells

**Recommendations** - Use extra wash step, handle plates carefully with lids on, use predetermined loading pattern before loading. Blot top of plate after rinsing.

### 3. *Color development very rapid; some color in healthy samples*

- a) Enzyme conjugate concentration too high
- b) Substrate concentration too high

**Recommendations** - Use enzyme conjugate and substrate concentrations that will give OD<sub>405 nm</sub> of about 1.0 in 30 to 60 min with good antigen source.