

1490: DESULFOVECTIS MEDIUM

NaCl	20.00	g
NH ₄ Cl	1.00	g
K ₂ HPO ₄	0.30	g
KH ₂ PO ₄	0.30	g
KCl	0.10	g
Na ₂ SO ₄	2.00	g
CaCl ₂ x 2 H ₂ O	0.10	g
Trace element solution SL-10	1.00	ml
Yeast extract	2.00	g
Sodium resazurin (0.1% w/v)	0.50	ml
Na ₂ CO ₃	1.50	g
Na ₂ S ₂ O ₃ x 5 H ₂ O	3.00	g
MgCl ₂ x 6 H ₂ O	3.00	g
L-Cysteine HCl x H ₂ O	0.50	g
Na ₂ S x 9 H ₂ O	0.50	g
Distilled water	1000.00	ml

1. Dissolve ingredients (except carbonate, thiosulfate, magnesium chloride, cysteine and sulfide), then sparge medium with 80% H₂ and 20% CO₂ gas mixture for 30 - 45 min to make it anoxic. Dispense medium under same gas atmosphere into anoxic Hungate-type tubes or serum vials to 30% of their volume and autoclave. Add magnesium chloride, cysteine and sulfide from sterile anoxic stock solutions prepared under 100% N₂ gas and carbonate from a sterile anoxic stock solution prepared under 80% N₂ and 20% CO₂ gas mixture. Adjust pH of the complete medium to 6.8 - 7.0, if necessary.
2. After inoculation pressurize cultivation vials to 1 bar overpressure with sterile 80% H₂ and 20% CO₂ gas mixture.

Trace element solution SL-10 (from medium 320)

HCl (25%)	10.00	ml
FeCl ₂ x 4 H ₂ O	1.50	g
ZnCl ₂	70.00	mg
MnCl ₂ x 4 H ₂ O	100.00	mg
H ₃ BO ₃	6.00	mg
CoCl ₂ x 6 H ₂ O	190.00	mg
CuCl ₂ x 2 H ₂ O	2.00	mg
NiCl ₂ x 6 H ₂ O	24.00	mg
Na ₂ MoO ₄ x 2 H ₂ O	36.00	mg
Distilled water	990.00	ml

First dissolve FeCl₂ in the HCl, then dilute in water, add and dissolve the other salts. Finally

make up to 1000.00 ml.