

836: DESULFITOBACTERIUM PCP-1 MEDIUM

KH ₂ PO ₄	0.27	g
K ₂ HPO ₄	0.35	g
NH ₄ Cl	0.53	g
MgCl ₂ x 6 H ₂ O	0.10	g
FeCl ₂ x 4 H ₂ O (0.1% w/v in 0.2 N HCl)	2.00	ml
CaCl ₂ x 2 H ₂ O	0.07	g
Trace element solution SL-10	1.00	ml
Yeast extract	1.00	g
Na-pyruvate	6.05	g
Sodium resazurin (0.1% w/v)	0.50	ml
Na ₂ CO ₃	0.50	g
Na ₂ SO ₃	0.63	g
Na ₂ S x 9 H ₂ O	0.05	g
Distilled water	1000.00	ml

Dissolve ingredients (except carbonate, sulfite and sulfide), adjust pH of solution to 7.0 and sparge with 100% N₂ gas for 30 - 45 min to make it anoxic. Dispense under same gas atmosphere into anoxic Hungate-type tubes or serum vials and autoclave. Prior to use add sulfite and sulfide from sterile anoxic stock solutions prepared under 100% N₂ gas and carbonate from a sterile anoxic stock solution prepared under 80% N₂ and 20% CO₂ gas atmosphere. The stock solution of sulfite should be freshly prepared and sterilized by filtration. Prior to use check pH of complete medium and adjust to 7.5, if necessary.

For DSM 13498: Supplement medium with 3.00 g/l Na-formate.

Trace element solution SL-10 (from medium 320)

HCl (25%)	10.00	ml
FeCl ₂ x 4 H ₂ O	1.50	g
ZnCl ₂	70.00	mg
MnCl ₂ x 4 H ₂ O	100.00	mg
H ₃ BO ₃	6.00	mg
CoCl ₂ x 6 H ₂ O	190.00	mg
CuCl ₂ x 2 H ₂ O	2.00	mg
NiCl ₂ x 6 H ₂ O	24.00	mg
Na ₂ MoO ₄ x 2 H ₂ O	36.00	mg
Distilled water	990.00	ml

First dissolve FeCl₂ in the HCl, then dilute in water, add and dissolve the other salts. Finally make up to 1000.00 ml.