

1096. UREAPLASMA MEDIUM**Solution A:**

PPLO-Broth	14.70	g
L-Cystein HCl	0.10	g
Fish-Sperm DNA	0.20	g
Phenol red	0.02	g
Urea	0.40	g
Water	700.00	ml

Solution B:

2% Yeastolate (BD)	100	ml
Horse serum	200	ml
Isovitalex	5	ml

Adjust solution A for final pH 6.0 with HCl. Autoclave solution A at 121°C for 15 min. Cool to room temperature and aseptically add solution B filter sterilised.