

111. BEIJERINCKIA MEDIUM

Glucose	10.0	g
K ₂ HPO ₄	0.8	g
KH ₂ PO ₄	0.2	g
MgSO ₄ x 7 H ₂ O	0.1	g
FeSO ₄ x 7 H ₂ O	20.0	mg
MnSO ₄ x 6 H ₂ O	2.0	mg
ZnSO ₄ x 6 H ₂ O	5.0	mg
CuSO ₄ x 6 H ₂ O	4.0	mg
Na ₂ MoO ₄ x 2 H ₂ O	5.0	mg
Agar	15.0	g
Distilled water	950.0	ml

Adjust pH to 6.5. Sterilize glucose separately (10 g in 50 ml H₂O) and mix after cooling.