

1113. LEPTOSPIRA-MEDIUM

Leptospira Medium Base EMJH (Difco)	2.3	g
agarose	1.5	g
water	900.0	ml
Leptospira Enrichment EMJH	100.0	ml

Dissolve the agarose in water by boiling. Add the Leptospira base and dissolve completely. Autoclave at 121°C for 15 min. Allow the medium to cool to 50-55°C. Warm the Leptospira Enrichment to approximately 25°C and add to the warm medium. Mix uniformly and dispense aseptically into sterile tubes as requested. Final pH is 7.5 at 25°C.