

1122. MEDIUM K

(NH ₄) ₂ SO ₄	2.000	g
KH ₂ PO ₄	2.000	g
NaCl	0.500	g
MgSO ₄ x 7 H ₂ O	0.125	g
FeSO ₄ x 7 H ₂ O	0.002	g
Agar Difco	20.000	g
Distilled water	1000.000	ml

Adjust to pH 7.2.

After autoclaving add methanol (sterilized separately) to a concentration of 1% (v/v).

For DSM 23053

+ 0.02 % Yeast extract (Difco)

Add 0.5% Methanol after autoclaving.

For DSM 25212, for DSM 25216 and for DSM 25045

Add 2.0% Methanol after autoclaving.

For DSM 23053

Add 0.02% Yeast extract.

For DSM 103579

Add 0.5% Methanol after autoclaving.