

1133. ENRICHED CYTOPHAGA AGAR MEDIUM

Tryptone	2.0	g
Beef extract	0.5	g
Yeast extract	0.5	g
Sodium acetate	0.2	g
Agar	15.0	g
Distilled water	1000.0	ml

Adjust to pH 7.2-7.4 and autoclave at 121°C/15 min.

For soft agar for the maintenance of active cultures reduce agar content to 4.0 g/l.
Dispense in 3 to 4 ml amounts in 7 ml screw-capped bottles (bijou bottles).