

**1141. THERMOGYMNOMONAS MEDIUM**

(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	0.20	g
KH <sub>2</sub> PO <sub>4</sub>	3.00	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	0.50	g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	0.25	g
Yeast extract (Difco)	1.00	g
Glucose	10.00	g
Distilled water	1000.00	ml

Adjust pH to 3.0 with 10 N H<sub>2</sub>SO<sub>4</sub>. Separately autoclave yeast extract and glucose as 10% solutions, and add aseptically to the medium.