

## 1157. PYEM MEDIUM

Peptone	2.0	g
Yeast extract	2.0	g
NH <sub>4</sub> Cl	0.5	g
Distilled water	1000.0	ml

After autoclaving and cooling add:

Riboflavin (0.2 mg/ml) filter-sterilized	5.0	ml
50% Glucose (autoclaved)	2.0	ml
20% MgSO <sub>4</sub> (autoclaved)	1.0	ml
10% CaCl <sub>2</sub> (autoclaved)	1.0	ml

Final pH 7.0