

1188. SPIROCHAETA CELLOBIOSIPHILA MEDIUM

Sea Salts (SIGMA)	35.0	g
Wolfe's Mineral Elixir (see medium 792)	1.0	ml
Yeast extract	1.0	g
Trypticase peptone (BD BBL)	0.5	g
Cellobiose	2.0	g
Distilled water	1000.0	ml

Dissolve compounds except cellobiose, adjust pH to 7.5 and dispense medium under air atmosphere into Hungate-type tubes or serum vials to 30% volume and autoclave. Add cellobiose after autoclaving from a stock solution sterilized by filtration.