

1265. BACTEROIDES GALACTURONICUS MEDIUM

Sodium polygalacturonate	4.00	g
Trypticase	5.00	g
Yeast Extract	2.50	g
MgSO ₄ × 7 H ₂ O	2.50	g
CaCl ₂ × 2 H ₂ O	0.15	g
FeSO ₄ × 7 H ₂ O	20.00	mg
(NH ₄) ₂ SO ₄	1.40	g
L-Cysteine	1.00	g
Resazurin	1.00	mg
NaHCO ₃	2.00	g
Distilled water	1000.00	ml

Adjust pH to 7.1 and prepare medium in a gas atmosphere of 100% N₂.

NaHCO₃ is prepared as 5% (w/v) solution, filter sterilized and added to the autoclaved medium.