

13. BACILLUS ACIDOCALDARIUS MEDIUM

Solution A:

Yeast extract	1.00	g
(NH ₄) ₂ SO ₄	0.20	g
MgSO ₄ x 7 H ₂ O	0.50	g
CaCl ₂ x 2 H ₂ O	0.25	g
KH ₂ PO ₄	0.60	g
Distilled water	500.00	ml
Adjust pH to 3.0 - 4.0.		

Solution B:

Glucose	1.00	g
Agar	20.00	g
Distilled water	500.00	ml

Sterilize solution A and B separately at 121°C for 15 min. Cool to 50°C and then combine.