

## 1335. MODIFIED T5 MEDIUM

Tryptone	0.5	g
Yeast extract	2.0	g
Glucose	1.0	g
Soluble starch	1.0	g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	1.0	g
Distilled water	1000.0	ml

Adjust to pH 7.0-7.5 and sterilize by autoclaving. The medium may be solidified by adding 20.0 g/l agar.

### For DSM 27250

Adjust to pH 7.5  $\pm$ 0.2