

1349. TG 408 MEDIUM

Mineral Salts Solution:

Peptone	5.0	g
Yeast Extract	0.5	g
Na-Pyruvate	10.0	g
Vitamin Solution 141	10.0	ml
Artificial Sea Water	750.0	ml
Distilled Water	250.0	ml

Vitamine Solution 141:

Biotin	2.0	mg
Folic acid	2.0	mg
Pyridoxine-HCl	10.0	mg
Thiamine-HCl x 2 H ₂ O	5.0	mg
Riboflavin	5.0	mg
Nicotinic acid	5.0	mg
D-Ca-pantothenate	5.0	mg
Vitamin B ₁₂	0.1	mg
p-Aminobenzoic acid	5.0	mg
Lipoic acid	5.0	mg
Distilled water	1000.0	ml

After autoclaving add 10 ml of Vitamine Solution 141 to 1 litre of the sterile Mineral Salts Solution. Final pH 7.5 - 8.0.