

1356. HMD MEDIUM**Mineral Salts Solution:**

Yeast Extract	0.10	g
Casamino acids	0.50	g
CaCl ₂ x 2 H ₂ O	0.10	g
MgCl ₂ x 6 H ₂ O	20.00	g
K ₂ SO ₄	5.00	g
NH ₄ Cl	0.50	g
KH ₂ PO ₄	0.05	g
NaCl	180.00	g
Distilled Water	1000.00	ml

Adjust to pH 7.5 and sterilize by autoclaving. The medium may be solidified by adding 12.0 g/l agar.