

1400. SALINIBACTER HALOPHILIC MEDIUM**Mineral Salts Solution:**

NaCl	200.0	g
MgSO ₄ x 7 H ₂ O	20.0	g
KH ₂ PO ₄	0.3	g
Yeast Extract	0.5	g
Bacto Peptone	0.5	g
Casamino Acids	0.5	g
Glucose	0.5	ml
Starch	0.5	g
Sodium Pyruvate	0.3	g
Distilled Water	1000.0	ml

Adjust to pH 7.0 and sterilize by autoclaving. The medium may be solidified by adding 20.0 g/l agar.