

**1429. SABOURAUD GLUCOSE MEDIUM**

SABOURAUD-2% Glucose-Bouillon (Merck 108339)	30.0	g
Agar	15.0	g
Distilled water	1000.0	ml

For DSM 6160, DSM 13536 after autoclaving add 3.2 ml/L of filter sterilized acetic acid (99% w/v).

The pH should be 3-4.

For DSM 5432 after autoclaving add 1.0 ml/L of filter sterilized acetic acid (99% w/v).

The pH should be 4-5.