

Glucose minerals salts medium

NaNO ₃	2.0	g
K ₂ HPO ₄	0.5	g
MgSO ₄ x 7 H ₂ O	0.2	g
MnSO ₄ x 5 H ₂ O	0.02	g
FeSO ₄ x 7 H ₂ O	0.02	g
CaCl ₂ x 7 H ₂ O	0.02	g
Glucose	5.0	g
Distilled water	1000.0	ml
Agar, if needed	15.0	g

Autoclave glucose separately. Adjust to pH 7.0