

1462: CAULOBACTER PYE MEDIUM

Peptone	2.00	g
Yeast extract (Difco)	1.00	g
CaCl ₂ x 2 H ₂ O	0.07	g
Agar, for solid medium (optional)	12.00	g
Distilled water	1000.00	ml

1. Add all the components together, except calcium chloride which is added from a sterilized stock solution when the autoclaved medium has cooled. Agar is added at a concentration of 12 g/l for solid media. For agar media, do not add supplemental calcium.
2. Adjust final pH to 7.0.