

1520. HALOPIGER MEDIUM (SG)

NaCl	250.0	g
KCl	3.0	g
Trisodium citrate	3.0	g
MgSO ₄ x 7 H ₂ O	20.0	g
Yeast extract (Difco)	1.0	g
Casamino acids (Difco)	7.5	g
Agar, if necessary	20.0	g

Add distilled water to give a final volume of 1000 ml.

Adjust pH to 8.0 before autoclaving!

If solid medium is prepared, autoclave agar and sodium chloride separately.