

## 1587. MSV MEDIUM FOR TETRACOCCUS CECHII

(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	0.500 g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	0.100 g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	0.050 g
K <sub>2</sub> HPO <sub>4</sub>	0.110 g
KH <sub>2</sub> PO <sub>4</sub>	0.085 g
FeCl <sub>3</sub> x 6 H <sub>2</sub> O	0.002 g
Na <sub>2</sub> -EDTA	0.003 g
Distilled water	1000.000 ml

After autoclaving the mineral base add the glucose at final concentration of 0.5g/l.

If solid media is required add agar to a final concentration of 15 g/l.