

172. CYTOPHAGA (marine) MEDIUM

Yeast extract (Difco)	1.0	g
Tryptone (Difco)	1.0	g
NaCl	24.7	g
KCl	0.7	g
MgSO ₄ x 7 H ₂ O	6.3	g
MgCl ₂ x 6 H ₂ O	4.6	g
CaCl ₂ x 2 H ₂ O	1.2	g
NaHCO ₃	0.2	g
Agar (Difco)	15.0	g
Distilled water	1000.0	ml

Adjust pH to 7.2. Sodium bicarbonate and calcium chloride are autoclaved separately, each in a small volume of distilled water.