

174. TREPONEMA STENOSTREPTUM MEDIUM

| | | |
|----------------------------------|--------|----|
| D-Glucose | 5.0 | g |
| Yeast extract (OXOID) | 2.0 | g |
| Peptone (BD BACTO) | 2.0 | g |
| Na-resazurin solution (0.1% w/v) | 0.5 | ml |
| Thioglycolate | 0.5 | g |
| Distilled water | 1000.0 | ml |

Dissolve ingredients, adjust pH to 7.6 with KOH and sparge medium with 100% N₂ gas for 30 – 45 min to make it anoxic. Distribute medium under same gas atmosphere into anoxic Hungate-type tubes or serum vials and autoclave.

For for stabs add 10.0 g/l agar.