

223. POLYANGIUM MEDIUM

KNO ₃	1.00	g
K ₂ HPO ₄	1.00	g
MgSO ₄ x 7 H ₂ O	0.30	g
CaCl ₂ x 2 H ₂ O	0.13	g
FeCl ₃ x 6 H ₂ O	32.00	mg
Agar	10.00	g
Tap water	1000.00	ml

Adjust pH to 7.2 (use pH indicator-paper).

Aseptically add a strip of filter paper to slants or broth. For plates, add 4 - 6 strips of filter paper to each surface.