

**226. AUREOBACTERIUM TERREGENS MEDIUM**

Casamino acids	2.0	g
Yeast extract	1.0	g
Glucose	1.0	g
(NH <sub>4</sub> ) citrate	1.0	g
K <sub>2</sub> HPO <sub>4</sub>	2.0	g
MgSO <sub>4</sub> × 7 H <sub>2</sub> O	0.5	g
FeCl <sub>3</sub> × 6 H <sub>2</sub> O	10.0	mg
Distilled water	1000.0	ml

Adjust pH to 7.0.

After sterilization, aseptically add 1 ml/l of a 10% solution of acetylacetone in ethanol.