

242. BMM - MEDIUM

$(\text{NH}_4)_2\text{SO}_4$	3.0	g
KH_2PO_4	2.0	g
K_2HPO_4	7.0	g
NaCl	2.0	g
$\text{MgSO}_4 \times 7 \text{ H}_2\text{O}$	0.5	g
Yeast extract	0.1	g
Biotin	10.0	μg
Thiamine-HCl $\times 2 \text{ H}_2\text{O}$	100.0	μg
$\text{MnSO}_4 \times \text{H}_2\text{O}$	6.2	mg
$\text{FeSO}_4 \times 7 \text{ H}_2\text{O}$	10.0	mg
Methanol	10.0	ml
Agar (if necessary)	15.0	g
Distilled water	1000.0	ml

Adjust pH to 7.2 prior to autoclaving.