

303. SIMONSIELLA MEDIUM

Trypticase (BBL)	17.0	g
Phytone (BBL)	3.0	g
NaCl	5.0	g
K ₂ HPO ₄	2.5	g
Yeast extract	4.0	g
Horse serum	100.0	ml
Agar, if necessary	15.0	g
Distilled water	1000.0	ml

Prepare the medium without horse serum. After autoclaving cool to 50°C, then add 100 ml/l of sterile horse serum. Transfer cultures when actively growing. Cultures die rapidly in stationary growth phase!