

**372. HALOBACTERIA MEDIUM**

Yeast extract	5.00	g
Casamino acids	5.00	g
Na-glutamate	1.00	g
KCl	2.00	g
Na <sub>3</sub> -citrate	3.00	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	20.00	g
NaCl	200.00	g
FeCl <sub>2</sub> x 4 H <sub>2</sub> O	36.00	mg
MnCl <sub>2</sub> x 4 H <sub>2</sub> O	0.36	mg
Agar	20.00	g

Add distilled water to give a final volume of 1000 ml.

Adjust pH to 7.0 - 7.2

For **DSM 102807** **DSM 102808** **DSM 102809** add 20g/l soluble starch.