

394. CHI 1776 MEDIUM

Tryptone	25.0	g
Yeast extract	7.5	g
1M Tris-HCl (pH 7.5)	20.0	ml
Distilled water	1000.0	ml

Autoclave, cool, and then add:

1M MgCl ₂	5.0	ml
1% Diaminopimelic acid	10.0	ml
0.4% Thymidine	10.0	ml
20% Glucose	25.0	ml

Magnesium chloride should be sterilized separately by autoclaving, and the other ingredients should be sterilized separately by filtration.