Microorganisms



395: DESULFUROCOCCUS MEDIUM

NH_4Cl 0.33	g
KH ₂ PO ₄ 0.33	g
KCI 0.33	g
$CaCl_2 \times 2 H_2O$ 0.44	g
MgCl2 x 6 H2O 0.70	g
NaCl 0.50	g
Trace element solution SL-10 1.00	ml
Yeast extract (OXOID) 0.20	g
Sodium resazurin (0.1% w/v) 0.50	ml
Sulfur (powdered) 10.00	g
D-Glucose 2.50	g
Wolin's vitamin solution (10x) 1.00	ml
$Na_2S \times 9 H_2O$ 0.50	g
Distilled water 1000.00	ml

Dissolve ingredients (except sulfur, glucose, vitamins and sulfide), then sparge medium with $80\%~N_2$ and $20\%~CO_2$ gas mixture for 30 - 45 min to make it anoxic. Adjust pH to 6.2 - 6.4 and dispense medium under $80\%~N_2$ and $20\%~CO_2$ gas atmosphere into anoxic Hungate-type tubes or serum vials containing already the appropriate amount of sulfur. Sterilize medium by heating cultivation vessels in a boiling water bath for 2 - 3 hours on each of 3 successive days. After sterilization add glucose, vitamins and sulfide from sterile anoxic stock solutions prepared under $100\%~N_2$ gas atmosphere. Vitamins are sterilized by filtration. Adjust pH of complete medium to 6.5, if necessary.

Trace element solution SL-10 (from medium 320)

HCI (25%)	10.00	ml
FeCl ₂ x 4 H ₂ O	1.50	g
ZnCl ₂	70.00	mg
$MnCl_2 \times 4 H_2O$	100.00	mg
H ₃ BO ₃	6.00	mg
CoCl ₂ x 6 H ₂ O	190.00	mg
CuCl ₂ x 2 H ₂ O	2.00	mg
$NiCl_2 \times 6 H_2O$	24.00	mg
$Na_2MoO_4 \times 2 H_2O$	36.00	mg
Distilled water	990.00	ml

First dissolve $FeCl_2$ in the HCl, then dilute in water, add and dissolve the other salts. Finally make up to 1000.00 ml.

Microorganisms

395: DESULFUROCOCCUS MEDIUM



Wolin's vitamin solution (10x) (from medium 120)

Biotin	20.00	mg
Folic acid	20.00	mg
Pyridoxine hydrochloride	100.00	mg
Thiamine HCI	50.00	mg
Riboflavin	50.00	mg
Nicotinic acid	50.00	mg
Calcium D-(+)-pantothenate	50.00	mg
Vitamin B ₁₂	1.00	mg
p-Aminobenzoic acid	50.00	mg
(DL)-alpha-Lipoic acid	50.00	mg
Distilled water	1000.00	ml