

430. ECTOTHIORHODOSPIRA ABDELMALEKII MEDIUM

KH ₂ PO ₄	0.80	g
NH ₄ Cl	0.80	g
CaCl ₂ x 2 H ₂ O	0.05	g
MgCl ₂ x 6 H ₂ O	0.10	g
NaCl	120.00	g
Na ₂ SO ₄	15.00	g
NaHCO ₃	10.00	g
Na ₂ CO ₃	4 - 5.00	g
Vitamin solution VA *	1.00	ml
Trace element solution SLA *	1.00	ml
Na ₂ S x 9 H ₂ O	0.50	g
Na-acetate	1.00	g
Distilled water	1000.00	ml

Adjust pH to 8.5, filter-sterilize and add appropriate amount of a sterile Na₂S solution.

* see medium 351