

431. ECTOTHIORHODOSPIRA VACUOLATA MEDIUM

KH ₂ PO ₄	1.00	g
NH ₄ Cl	0.50	g
MgCl ₂ x 6 H ₂ O	0.20	g
CaCl ₂ x 2 H ₂ O	0.10	g
NaHCO ₃	3.00	g
NaCl	30.00	g
Trace element solution SLA*	1.00	ml
Vitamin solution VA*	1.00	ml
Na ₂ S x 9 H ₂ O	0.05	g
Na ₂ S ₂ O ₃ x 5 H ₂ O	0.50	g
Na-malate	1.00	g
Distilled water	1000.00	ml

Adjust pH to 8.7, filter sterilize and add appropriate amount of Na₂S solution.

* see medium 351