

514b. MEDIUM 514 plus additional salt

Use medium 514 to which is added NaCl at the concentration indicated (i.e. "514 + 6% NaCl" indicates use medium 514 to which 60 g/l NaCl has been added). Final pH 7.6. The medium may be solidified by adding 15 - 20 g/l agar.

For strains DSM 21167, DSM 21194, DSM 23711 100.0 g

For strains DSM 5349, DSM 19800, DSM 24106, DSM 24107, DSM 24186, DSM 23100, DSM 21021, DSM 19800, DSM 19821, DSM 24314
NaCl 80.0 g/l

For strains DSM 23438, DSM 21056, DSM 21166, DSM 21888, DSM 5160, DSM 4683, DSM 4684, DSM 22872
NaCl 60.0 g/l

For strains DSM 5351, DSM 23229, DSM 22730, DSM 22404, DSM 22403, DSM 22402, DSM 25260
NaCl 50.0 g/l

For strains DSM 24014
NaCl 40.0 g/l

For strains DSM 23437, DSM 23508, DSM 23509, DSM 21076, DSM 26132
NaCl 30.0 g/l