

**659. M17 MEDIUM (modified)**

Casein peptone	5.00	g
Soya peptone	5.00	g
Meat peptone	5.00	g
Yeast extract	2.50	g
Ascorbic acid	0.50	g
MgSO <sub>4</sub> × 7 H <sub>2</sub> O	0.25	g
Na <sub>2</sub> -β-glycerolphosphate	9.50	g
Distilled water	1000.00	ml

After autoclaving add a lactose solution (sterilized by filtration) up to a concentration of 8.0 g/l, and 1.2 ml 1M CaCl<sub>2</sub> × 2 H<sub>2</sub>O. Adjust pH to 7.15 ± 0.05.