

**767. ERYTHROMICROBIUM AND ROSEOCOCCUS MEDIUM**

NH <sub>4</sub> Cl	0.30	g
K <sub>2</sub> HPO <sub>4</sub>	0.30	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	0.50	g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	0.05	g
KCl	0.30	g
Na-acetate	1.00	g
Yeast extract	1.00	g
Bactopeptone	1.00	g
Trace elements SL-6 (see medium 27)	1.00	ml
Vitamin B <sub>12</sub>	20.00	µg
Distilled water	1000.00	ml

Final pH 7.5-7.8

Add vitamin B<sub>12</sub> from a filter sterilised stock solution.

For **DSM 8509** add 15g /l NaCl.