

779. CAPNOCYTOPHAGA II MEDIUM

Lab-lemco meat extract (Oxoid)	2.4	g
Proteose peptone no.3 (Difco)	10.0	g
Yeast extract (Difco, Bacto)	5.0	g
Na ₂ HPO ₄	4.0	g
Glucose	1.5	g
Starch, soluble	0.5	g
Cysteine-HCl x H ₂ O	0.5	g
Agar, if necessary	15.0	g
Horse blood	50.0	ml
Distilled water	1000.0	ml

Adjust pH to 7.6 - 7.8.

Add the horse blood to the heat sterilized medium after being cooled to 50°C. Incubate under air plus 5% CO₂ or anaerobically under 95% N₂ plus 5% CO₂.