

805. METHYLOBACTERIUM THIOCYANATUM MEDIUM

KSCN	0.25	g
MgSO ₄ x 7 H ₂ O	0.10	g
K ₂ HPO ₄	1.50	g
Na ₂ HPO ₄ x 2 H ₂ O	7.90	g
Glucose	4.50	g
2% FeSO ₄ x 7 H ₂ O	1.00	ml
Distilled water	1000.00	ml

To solid media add 15 g/l agar. Sterilise medium at 115°C for 10 minutes.