

84. ROLLED OATS MINERAL MEDIUM

Rolled oats	20.0	g
Agar	12.0	g
Trace element solution (see below)	1.0	ml
Distilled water	1000.0	ml

Rolled oats are boiled for 20 min and filtered, then add agar. Make up to 1000 ml. Then trace elements are added:

Trace element solution:

FeSO ₄ x 7 H ₂ O	0.1	g
MnCl ₂ x 4 H ₂ O	0.1	g
ZnSO ₄ x 7 H ₂ O	0.1	g
Distilled water	100.0	ml