

## 852. GYP GLUCOSE-YEAST-PEPTONE MEDIUM

Glucose	20	g
Yeast extract	10	g
Peptone (tryptic)	10	g
Na-acetate	10	g
Salt solution (see below)	5	ml
Distilled water	1000	ml

*Salt solution:*

MgSO <sub>4</sub> x 7 H <sub>2</sub> O	40	g
MnSO <sub>4</sub> x 4 H <sub>2</sub> O	2	g
FeSO <sub>4</sub> x 7 H <sub>2</sub> O	2	g
NaCl	2	g
Distilled water	1000	ml

For solid media add 10 g/l bacteriological agar. Adjust to pH 6.8. After sterilization of the medium add filter sterilized salt solution.