

87a. CHLOROFLEXUS AGGREGANS MEDIUM

To medium 87 add 1.0 ml/litre sterile vitamin solution CA.

Vitamin solution CA

Distilled water	100.0	ml
Nicotinic acid	100.0	mg
Thiamine-HCl x 2 H ₂ O	100.0	mg
Biotin	5.0	mg
p-Aminobenzoic acid	50.0	mg
Vitamin B ₁₂	1.0	mg
Ca-pantothenate	50.0	mg
Pyridoxine-HCl	50.0	mg
Folic acid	50.0	mg
Na ₃ -EDTA	200.0	mg

Adjust pH to 7.5. The solution is filter-sterilized.