

888. AQUASPIRILLUM MEDIUM

(NH ₄) ₂ SO ₄	1.0	g
MgSO ₄ x 7 H ₂ O	1.0	g
CaCl ₂ x 6 H ₂ O	30.0	mg
Na ₂ H ₂ PO ₄	10.0	mg
Sodium succinate	1.00	g
Casamino acids (DIFCO)	1.5	g
Trace element solution (see medium 1369)	1.0	ml
Agar (DIFCO)	0.5	g
Distilled water	1000.0	ml
pH 7.5		

Add after autoclaving sterile solutions of

Na ₂ S ₂ O ₃ x 5 H ₂ O (10% solution)	1.0	ml
vitamin solution (see medium 1478)	1.0	ml