

910. MODIFIED CM+YE MEDIUM (B)

Bacto vitamin assay		
casamino acids (Difco)	7.50	g
Bacto yeast extract (Difco)	10.00	g
MgSO ₄ x 7 H ₂ O	10.00	g
Na ₃ Citrate x 2 H ₂ O	3.00	g
KCl	2.00	g
NaCl	100.00	g
Fe ²⁺ solution*	1.00	ml
Agar	15.00	g
Distilled water	1000.00	ml

Adjust pH to 7.4

*Fe²⁺ solution:

FeSO ₄ x 7 H ₂ O	4.98	g
Distilled water	100.00	ml