

**966. R3 MEDIUM**

Yeast extract	1.00	g
Proteose peptone (Difco No.3)	1.00	g
Casamino acids (Difco)	1.00	g
Glucose	1.00	g
K <sub>2</sub> HPO <sub>4</sub>	0.60	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	0.10	g
Na pyruvate	0.05	g
Distilled water	1000.00	ml

Final pH should be 7.8

For **DSM 17972** prepare medium without glucose and adjust to pH 8.0 with 1M NaOH. Autoclave at 121°C for 15 minutes.