

97. HALOBACTERIUM MEDIUM

Casamino acids	7.50	g
Yeast extract	10.00	g
Na ₃ -citrate	3.00	g
KCl	2.00	g
MgSO ₄ x 7 H ₂ O	20.00	g
FeSO ₄ x 7 H ₂ O	0.05	g
MnSO ₄ x H ₂ O	0.20	mg
NaCl	250.00	g
Agar	20.00	g
Distilled water	1000.00	ml

Adjust pH to 7.4. Add the agar after dissolving all ingredients in the water and adjustment of pH.