

**1122a. METHYLOVERSATILIS THERMOTOLERANS MEDIUM**

(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	2.000	g
KH <sub>2</sub> PO <sub>4</sub>	2.000	g
NaCl	0.500	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	0.025	g
FeSO <sub>4</sub> x 7 H <sub>2</sub> O	0.002	g
Agar Difco	20.000	g
Distilled water	1000.000	ml

Adjust to pH 7.2.

After autoclaving add methanol (sterilized separately) to a concentration of 0.5 % (v/v) and 0.1 ml vitamin solution B<sub>12</sub> (10µg /ml, filter sterilized).