

**1220. ACIDOSOMA MEDIUM**

KH <sub>2</sub> PO <sub>4</sub>	0.10	g
(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	0.25	g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	0.05	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	0.02	g
Trace elements (see medium 922)	1.00	ml
Yeast extract	0.10	g
Na gluconate	0.50	g
Distilled water	1000.00	ml

Adjust to pH 5.0 – 5.5. Agar may be added at 15 g/l for solid media.